



KLRAP1 blocking peptide (CDBP5647)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	KLRAP1 killer cell lectin-like receptor subfamily A pseudogene 1 [Homo sapiens (human)]
Official Symbol	KLRAP1
Synonyms	KLRAP1; killer cell lectin-like receptor subfamily A pseudogene 1; Ly49; KLRA1; LY49L; Ly-49L
Entrez Gene ID	10748