



## Klra3 blocking peptide (CDBP5650)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

### GENE INFORMATION

|                 |                                                                                                                                                                                                                                                |
|-----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">Klra3 killer cell lectin-like receptor, subfamily A, member 3 [ Mus musculus (house mouse) ]</a>                                                                                                                                   |
| Official Symbol | Klra3                                                                                                                                                                                                                                          |
| Synonyms        | Klra3; killer cell lectin-like receptor, subfamily A, member 3; 5E6; Nk2; Nk-2; Ly49c; Nk2.1; NK-2.1; killer cell lectin-like receptor 3; lymphocyte antigen 49c; T-cell surface glycoprotein Ly-49C; natural killer cell-associated antigen 2 |
| Entrez Gene ID  | <a href="#">16634</a>                                                                                                                                                                                                                          |
| mRNA Refseq     | <a href="#">NM_001289604</a>                                                                                                                                                                                                                   |
| Protein Refseq  | <a href="#">NP_001276533</a>                                                                                                                                                                                                                   |
| UniProt ID      | Q64329                                                                                                                                                                                                                                         |

