



KPNA7 blocking peptide (CDBP5656)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	KPNA7 karyopherin alpha 7 (importin alpha 8) [Homo sapiens (human)]
Official Symbol	KPNA7
Synonyms	KPNA7; karyopherin alpha 7 (importin alpha 8); importin subunit alpha-8; karyopherin 7; karyopherin subunit alpha-7
Entrez Gene ID	402569
mRNA Refseq	NM_001145715
Protein Refseq	NP_001139187
UniProt ID	A9QM74
Pathway	Antiviral mechanism by IFN-stimulated genes; Cytokine Signaling in Immune system; ISG15

antiviral mechanism; Immune System; Interferon Signaling

Function

protein transporter activity
