



L1TD1 blocking peptide (CDBP5659)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	L1TD1 LINE-1 type transposase domain containing 1 [Homo sapiens (human)]
Official Symbol	L1TD1
Synonyms	L1TD1; LINE-1 type transposase domain containing 1; ECAT11; LINE-1 type transposase domain-containing protein 1; ES cell-associated protein 11
Entrez Gene ID	54596
mRNA Refseq	NM_001164835
Protein Refseq	NP_001158307
UniProt ID	Q5T7N2