



Lamp2 blocking peptide (CDBP5662)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	50 µg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Lamp2 lysosomal-associated membrane protein 2 [Mus musculus (house mouse)]
Official Symbol	Lamp2
Synonyms	Lamp2; lysosomal-associated membrane protein 2; Mac3; LGP-B; CD107b; Lamp-2; Lamp II; Lamp-2a; Lamp-2b; Lamp-2c; lysosome-associated membrane glycoprotein 2; lysosomal membrane glycoprotein 2; CD107 antigen-like family member B; lysosomal membrane glycoprotein type B
Entrez Gene ID	16784
mRNA Refseq	NM_001017959
Protein Refseq	NP_001017959
UniProt ID	P17047

Pathway	Hemostasis; Lysosome; Phagosome; Platelet activation; Platelet degranulation; Response to elevated platelet cytosolic Ca ²⁺ ; Senescence and Autophagy; Tuberculosis
Function	enzyme binding; protein domain specific binding