



LAMTOR1 blocking peptide (CDBP5663)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LAMTOR1 late endosomal/lysosomal adaptor, MAPK and MTOR activator 1 [Homo sapiens (human)]
Official Symbol	LAMTOR1
Synonyms	LAMTOR1; late endosomal/lysosomal adaptor, MAPK and MTOR activator 1; p18; PDRO; C11orf59; p27RF-Rho; Ragulator1; ragulator complex protein LAMTOR1; RhoA activator C11orf59; lipid raft adaptor protein p18; ragulator complex protein PDRO; p27Kip1-releasing factor from RhoA; p27kip1 releasing factor from RhoA; protein associated with DRMs and endosomes; late endosomal/lysosomal adaptor and MAPK and MTOR activator 1
Entrez Gene ID	55004
mRNA Refseq	NM_017907
Protein Refseq	NP_060377

UniProt ID	Q6IAA8
Function	contributes_to guanyl-nucleotide exchange factor activity; protein binding; contributes_to protein complex scaffold