



LRRC1 blocking peptide (CDBP5664)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRRC1 leucine rich repeat containing 1 [Homo sapiens (human)]
Official Symbol	LRRC1
Synonyms	LRRC1; leucine rich repeat containing 1; LANO; dJ523E19.1; leucine-rich repeat-containing protein 1; LANO adapter protein; LAP and no PDZ protein; leucine-rich repeat-containing 1; LAP (leucine-rich repeats and PDZ) and no PDZ protein
Entrez Gene ID	55227
mRNA Refseq	NM_018214
Protein Refseq	NP_060684
UniProt ID	Q9BTT6