



CERS6 blocking peptide (CDBP5666)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	CERS6 ceramide synthase 6 [Homo sapiens (human)]
Official Symbol	CERS6
Synonyms	CERS6; ceramide synthase 6; CERS5; LASS6; longevity assurance homolog 6; LAG1 homolog, ceramide synthase 6
Entrez Gene ID	253782
mRNA Refseq	NM_001256126
Protein Refseq	NP_001243055
UniProt ID	Q6ZMG9
Pathway	Ceramide biosynthesis; Metabolism; Metabolism of lipids and lipoproteins; Sphingolipid

Metabolism; Sphingolipid de novo biosynthesis; Sphingolipid metabolism; Sphingosine biosynthesis

Function

sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity; sphingosine N-acyltransferase activity
