



LEMD1 blocking peptide (CDBP5669)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LEMD1 LEM domain containing 1 [Homo sapiens (human)]
Official Symbol	LEMD1
Synonyms	LEMD1; LEM domain containing 1; CT50; LEMP-1; LEM domain-containing protein 1; LEM domain protein 1; cancer/testis antigen 50
Entrez Gene ID	93273
mRNA Refseq	NM_001001552
Protein Refseq	NP_001001552
UniProt ID	Q68G75