



LEMD2 blocking peptide (CDBP5670)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LEMD2 LEM domain containing 2 [Homo sapiens (human)]
Official Symbol	LEMD2
Synonyms	LEMD2; LEM domain containing 2; NET25; dJ482C21.1; LEM domain-containing protein 2; hLEM2
Entrez Gene ID	221496
mRNA Refseq	NM_001143944
Protein Refseq	NP_001137416
UniProt ID	Q8NC56
Pathway	Cell Cycle; Clearance of Nuclear Envelope Membranes from Chromatin; Depolymerisation of

the Nuclear Lamina; Initiation of Nuclear Envelope Reformation; M Phase; Mitotic Anaphase;
Mitotic Metaphase and Anaphase; Mitotic Prophase
