



MORF4 blocking peptide (CDBP5671)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	FAIM2 Fas apoptotic inhibitory molecule 2 [Homo sapiens (human)]
Official Symbol	MORF4
Synonyms	FAIM2; Fas apoptotic inhibitory molecule 2; LFG; LFG2; NGP35; NMP35; TMBIM2; protein lifeguard 2; neural membrane protein 35; transmembrane BAX inhibitor motif-containing protein 2
Entrez Gene ID	23017
mRNA Refseq	NM_012306
Protein Refseq	NP_036438
UniProt ID	Q9BWQ8

