



ST8SIA4 blocking peptide (CDBP5675)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	Lgi2 leucine-rich repeat LGI family, member 2 [Mus musculus (house mouse)]
Official Symbol	ST8SIA4
Synonyms	Lgi2; leucine-rich repeat LGI family, member 2; mKIAA1916; leucine-rich repeat LGI family member 2; leucine-rich glioma-inactivated protein 2
Entrez Gene ID	246316
mRNA Refseq	NM_144945
Protein Refseq	NP_659194
UniProt ID	Q8K4Z0
Function	molecular_function