



PTPRD blocking peptide (CDBP5676)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LGI3 leucine-rich repeat LGI family, member 3 [Homo sapiens (human)]
Official Symbol	PTPRD
Synonyms	LGI3; leucine-rich repeat LGI family, member 3; LGIL4; leucine-rich repeat LGI family member 3; LGI1-like protein 4; leucine-rich glioma-inactivated protein 3
Entrez Gene ID	203190
mRNA Refseq	NM_139278
Protein Refseq	NP_644807
UniProt ID	Q8N145
Function	catalytic activity