



LGI4 blocking peptide (CDBP5677)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LGI4 leucine-rich repeat LGI family, member 4 [Homo sapiens (human)]
Official Symbol	LGI4
Synonyms	LGI4; leucine-rich repeat LGI family, member 4; LGIL3; leucine-rich repeat LGI family member 4; LGI1-like protein 3; leucine-rich glioma-inactivated gene 4; leucine-rich glioma-inactivated protein 4
Entrez Gene ID	163175
mRNA Refseq	NM_139284
Protein Refseq	NP_644813
UniProt ID	Q8N135