



## LGI4 blocking peptide (CDBP5678)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

### GENE INFORMATION

Gene Name	<a href="#">LGI4 leucine-rich repeat LGI family, member 4 [ Homo sapiens (human) ]</a>
Official Symbol	LGI4
Synonyms	LGI4; leucine-rich repeat LGI family, member 4; LGIL3; leucine-rich repeat LGI family member 4; LGI1-like protein 3; leucine-rich glioma-inactivated gene 4; leucine-rich glioma-inactivated protein 4
Entrez Gene ID	<a href="#">163175</a>
mRNA Refseq	<a href="#">NM_139284</a>
Protein Refseq	<a href="#">NP_644813</a>
UniProt ID	Q8N135