



ANKRD54 blocking peptide (CDBP5680)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ANKRD54 ankyrin repeat domain 54 [Homo sapiens (human)]
Official Symbol	ANKRD54
Synonyms	ANKRD54; ankyrin repeat domain 54; LIAR; ankyrin repeat domain-containing protein 54; lyn-interacting ankyrin repeat protein
Entrez Gene ID	129138
mRNA Refseq	NM_138797
Protein Refseq	NP_620152
UniProt ID	Q6NXT1
Function	protein complex binding; protein kinase regulator activity