



LRFN1 blocking peptide (CDBP5690)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRFN1 leucine rich repeat and fibronectin type III domain containing 1 [Homo sapiens (human)]
Official Symbol	LRFN1
Synonyms	LRFN1; leucine rich repeat and fibronectin type III domain containing 1; SALM2; leucine-rich repeat and fibronectin type III domain-containing protein 1; synaptic adhesion-like molecule 2
Entrez Gene ID	57622
mRNA Refseq	NM_020862
Protein Refseq	NP_065913
UniProt ID	Q9P244