



LRFN2 blocking peptide (CDBP5691)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRFN2 leucine rich repeat and fibronectin type III domain containing 2 [Homo sapiens (human)]
Official Symbol	LRFN2
Synonyms	LRFN2; leucine rich repeat and fibronectin type III domain containing 2; SALM1; FIGLER2; KIAA1246; leucine-rich repeat and fibronectin type-III domain-containing protein 2; synaptic adhesion-like molecule 1; fibronectin type III, immunoglobulin and leucine rich repeat domains 2
Entrez Gene ID	57497
mRNA Refseq	NM_020737
Protein Refseq	NP_065788
UniProt ID	Q9ULH4