



LRFN3 blocking peptide (CDBP5692)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRFN3 leucine rich repeat and fibronectin type III domain containing 3 [Homo sapiens (human)]
Official Symbol	LRFN3
Synonyms	LRFN3; leucine rich repeat and fibronectin type III domain containing 3; SALM4; FIGLER1; leucine-rich repeat and fibronectin type-III domain-containing protein 3; synaptic adhesion-like molecule 4; fibronectin type III, immunoglobulin and leucine rich repeat domains 1
Entrez Gene ID	79414
mRNA Refseq	NM_024509
Protein Refseq	NP_078785
UniProt ID	Q9BTN0