



LRFN4 blocking peptide (CDBP5693)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRFN4 leucine rich repeat and fibronectin type III domain containing 4 [Homo sapiens (human)]
Official Symbol	LRFN4
Synonyms	LRFN4; leucine rich repeat and fibronectin type III domain containing 4; SALM3; SALM3.; FIGLER6; leucine-rich repeat and fibronectin type-III domain-containing protein 4; fibronectin type III, immunoglobulin and leucine rich repeat domains 6
Entrez Gene ID	78999
mRNA Refseq	NM_024036
Protein Refseq	NP_076941
UniProt ID	Q6PJG9