



LRFN5 blocking peptide (CDBP5694)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRFN5 leucine rich repeat and fibronectin type III domain containing 5 [Homo sapiens (human)]
Official Symbol	LRFN5
Synonyms	LRFN5; leucine rich repeat and fibronectin type III domain containing 5; SALM5; FIGLER8; C14orf146; leucine-rich repeat and fibronectin type-III domain-containing protein 5; fibronectin type III, immunoglobulin and leucine rich repeat domains 8
Entrez Gene ID	145581
mRNA Refseq	NM_152447
Protein Refseq	NP_689660
UniProt ID	Q96NI6