



LRRC19 blocking peptide (CDBP5695)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRRC19 leucine rich repeat containing 19 [Homo sapiens (human)]
Official Symbol	LRRC19
Synonyms	LRRC19; leucine rich repeat containing 19; leucine-rich repeat-containing protein 19
Entrez Gene ID	64922
mRNA Refseq	NM_022901
Protein Refseq	NP_075052
UniProt ID	Q9H756