



LRRC26 blocking peptide (CDBP5696)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRRC26 leucine rich repeat containing 26 [Homo sapiens (human)]
Official Symbol	LRRC26
Synonyms	LRRC26; leucine rich repeat containing 26; CAPC; bA350O14.10; leucine-rich repeat-containing protein 26; cytokeratin associated protein; cytokeratin-associated protein in cancer; BK channel auxiliary gamma subunit LRRC26; BK channel auxiliary gamma subunit LRRC26
Entrez Gene ID	389816
mRNA Refseq	NM_001013653
Protein Refseq	NP_001013675
UniProt ID	Q2I0M4

Function

potassium channel regulator activity; protein binding
