



LRRFIP2 blocking peptide (CDBP5697)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRRFIP2 leucine rich repeat (in FLII) interacting protein 2 [Homo sapiens (human)]
Official Symbol	LRRFIP2
Synonyms	LRRFIP2; leucine rich repeat (in FLII) interacting protein 2; HUFI-2; leucine-rich repeat flightless-interacting protein 2; LRR FLII-interacting protein 2
Entrez Gene ID	9209
mRNA Refseq	NM_001134369
Protein Refseq	NP_001127841
UniProt ID	Q9Y608
Function	LRR domain binding; protein binding