



LRRTM1 blocking peptide (CDBP5698)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRRTM1 leucine rich repeat transmembrane neuronal 1 [Homo sapiens (human)]
Official Symbol	LRRTM1
Synonyms	LRRTM1; leucine rich repeat transmembrane neuronal 1; leucine-rich repeat transmembrane neuronal protein 1; leucine-rich repeat transmembrane neuronal 1 protein
Entrez Gene ID	347730
mRNA Refseq	NM_178839
Protein Refseq	NP_849161
UniProt ID	Q86UE6
Function	molecular_function