



LRRTM2 blocking peptide (CDBP5699)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRRTM2 leucine rich repeat transmembrane neuronal 2 [Homo sapiens (human)]
Official Symbol	LRRTM2
Synonyms	LRRTM2; leucine rich repeat transmembrane neuronal 2; leucine-rich repeat transmembrane neuronal protein 2; leucine-rich repeat neuronal 2 protein; leucine-rich repeat transmembrane neuronal 2 protein
Entrez Gene ID	26045
mRNA Refseq	NM_015564
Protein Refseq	NP_056379
UniProt ID	O43300

Functionneurexin family protein binding
