



LRRTM3 blocking peptide (CDBP5700)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LRRTM3 leucine rich repeat transmembrane neuronal 3 [Homo sapiens (human)]
Official Symbol	LRRTM3
Synonyms	LRRTM3; leucine rich repeat transmembrane neuronal 3; leucine-rich repeat transmembrane neuronal protein 3
Entrez Gene ID	347731
mRNA Refseq	NM_178011
Protein Refseq	NP_821079
UniProt ID	Q86VH5