



LYRM2 blocking peptide (CDBP5708)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LYRM2 LYR motif containing 2 [Homo sapiens (human)]
Official Symbol	LYRM2
Synonyms	LYRM2; LYR motif containing 2; DJ122O8.2; LYR motif-containing protein 2
Entrez Gene ID	57226
mRNA Refseq	NM_020466
Protein Refseq	NP_065199
UniProt ID	Q9NU23