



NDUFB9 blocking peptide (CDBP5709)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NDUFB9 NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 9, 22kDa [Homo sapiens (human)]
Official Symbol	NDUFB9
Synonyms	NDUFB9; NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 9, 22kDa; B22; LYRM3; CI-B22; UQOR22; NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 9; complex I-B22; complex I B22 subunit; LYR motif-containing protein 3; NADH-ubiquinone oxidoreductase B22 subunit
Entrez Gene ID	4715
mRNA Refseq	NM_001278645
Protein Refseq	NP_001265574

UniProt ID	Q9Y6M9
Pathway	Alzheimer's disease; Electron Transport Chain; Huntington's disease; Metabolism; NADH dehydrogenase (ubiquinone) 1 beta subcomplex; Non-alcoholic fatty liver disease (NAFLD); Oxidative phosphorylation; Parkinson's disease
Function	NADH dehydrogenase (ubiquinone) activity