



LYZL1 blocking peptide (CDBP5710)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LYZL1 lysozyme-like 1 [Homo sapiens (human)]
Official Symbol	LYZL1
Synonyms	LYZL1; lysozyme-like 1; LYC2; KAAG648; PRO1278; bA534G20.1; lysozyme-like protein 1
Entrez Gene ID	84569
mRNA Refseq	NM_032517
Protein Refseq	NP_115906
UniProt ID	Q6UWQ5
Function	lysozyme activity