



SEC16B blocking peptide (CDBP5712)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SEC16B is a mammalian homolog of <i>S. cerevisiae</i> Sec16 that is required for organization of transitional endoplasmic reticulum (ER) sites and protein export (Bhattacharyya and Glick, 2007 [PubMed 17192411]).[supplied by OMIM, Jun 2009]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	SEC16B SEC16 homolog B (S. cerevisiae) [Homo sapiens (human)]
Official Symbol	SEC16B
Synonyms	SEC16B; SEC16 homolog B (<i>S. cerevisiae</i>); RGPR; LZTR2; SEC16S; PGPR-p117; protein transport protein Sec16B; RGPR-p117; protein SEC16 homolog B; leucine zipper transcription regulator 2; regucalcin gene promotor region related protein; regucalcin gene promoter region-related protein p117
Entrez Gene ID	89866
mRNA Refseq	NM_033127

Protein Refseq [NP_149118](#)

UniProt ID Q96JE7
