



LZTS2 blocking peptide (CDBP5714)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LZTS2 leucine zipper, putative tumor suppressor 2 [Homo sapiens (human)]
Official Symbol	LZTS2
Synonyms	LZTS2; leucine zipper, putative tumor suppressor 2; LAPSER1; leucine zipper putative tumor suppressor 2; hLZTS2
Entrez Gene ID	84445
mRNA Refseq	NM_032429
Protein Refseq	NP_115805
UniProt ID	Q9BRK4
Function	protein binding