



NAA35 blocking peptide (CDBP5719)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NAA35 N(alpha)-acetyltransferase 35, NatC auxiliary subunit [Homo sapiens (human)]
Official Symbol	NAA35
Synonyms	NAA35; N(alpha)-acetyltransferase 35, NatC auxiliary subunit; EGAP; MAK10; MAK10P; bA379P1.1; N-alpha-acetyltransferase 35, NatC auxiliary subunit; protein MAK10 homolog; corneal wound healing-related protein; embryonic growth-associated protein homolog; MAK10 homolog, amino-acid N-acetyltransferase subunit
Entrez Gene ID	60560
mRNA Refseq	NM_024635
Protein Refseq	NP_078911
UniProt ID	Q5VZE5

Functionprotein binding
