



MATN2 blocking peptide (CDBP5731)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a member of the von Willebrand factor A domain containing protein family. This family of proteins is thought to be involved in the formation of filamentous networks in the extracellular matrices of various tissues. This protein contains five von Willebrand factor A domains. The specific function of this gene has not yet been determined. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MATN2 matrilin 2 [Homo sapiens (human)]
Official Symbol	MATN2
Synonyms	MATN2; matrilin 2; matrilin-2
Entrez Gene ID	4147
mRNA Refseq	NM_002380

Protein Refseq	NP_002371
UniProt ID	O00339
Function	calcium ion binding; molecular_function; protein binding