



LY86 blocking peptide (CDBP5741)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	LY86 lymphocyte antigen 86 [Homo sapiens (human)]
Official Symbol	LY86
Synonyms	LY86; lymphocyte antigen 86; MD-1; MMD-1; dJ80N2.1; ly-86; protein MD-1; MD-1, RP105-associated
Entrez Gene ID	9450
mRNA Refseq	NM_004271
Protein Refseq	NP_004262
UniProt ID	O95711
Pathway	Immune System; Innate Immune System; Toll Like Receptor 4 (TLR4) Cascade; Toll-Like

