



MED28 blocking peptide (CDBP5742)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MED28 mediator complex subunit 28 [Homo sapiens (human)]
Official Symbol	MED28
Synonyms	MED28; mediator complex subunit 28; EG1; magicin; 1500003D12Rik; mediator of RNA polymerase II transcription subunit 28; endothelial-derived gene 1; endothelial-derived protein 1; tumor angiogenesis marker EG-1; merlin and Grb2-interacting cytoskeletal protein; mediator of RNA polymerase II transcription, subunit 28 homolog
Entrez Gene ID	80306
mRNA Refseq	NM_025205
Protein Refseq	NP_079481
UniProt ID	Q9H204

Functionactin binding; protein binding
