



MEIG1 blocking peptide (CDBP5744)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MEIG1 meiosis/spermiogenesis associated 1 [Homo sapiens (human)]
Official Symbol	MEIG1
Synonyms	MEIG1; meiosis/spermiogenesis associated 1; SPATA39; bA2K17.3; meiosis expressed gene 1 protein homolog; spermatogenesis associated 39; meiosis expressed gene 1 homolog
Entrez Gene ID	644890
mRNA Refseq	NM_001080836
Protein Refseq	NP_001074305
UniProt ID	Q5JSS6