



METTL7A blocking peptide (CDBP5745)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	METTL7A methyltransferase like 7A [Homo sapiens (human)]
Official Symbol	METTL7A
Synonyms	METTL7A; methyltransferase like 7A; AAM-B; methyltransferase-like protein 7A; protein AAM-B
Entrez Gene ID	25840
mRNA Refseq	NM_014033
Protein Refseq	NP_054752
UniProt ID	Q9H8H3
Function	methyltransferase activity