



METTL7B blocking peptide (CDBP5748)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	METTL7B methyltransferase like 7B [Homo sapiens (human)]
Official Symbol	METTL7B
Synonyms	METTL7B; methyltransferase like 7B; ALDI; methyltransferase-like protein 7B; associated with lipid droplets 1
Entrez Gene ID	196410
mRNA Refseq	NM_152637
Protein Refseq	NP_689850
UniProt ID	Q6UX53
Function	methyltransferase activity