



MFSD1 blocking peptide (CDBP5750)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MFSD1 major facilitator superfamily domain containing 1 [Homo sapiens (human)]
Official Symbol	MFSD1
Synonyms	MFSD1; major facilitator superfamily domain containing 1; SMAP4; major facilitator superfamily domain-containing protein 1; smooth muscle cell-associated protein 4
Entrez Gene ID	64747
mRNA Refseq	NM_001167903
Protein Refseq	NP_001161375
UniProt ID	Q9H3U5