



mfsd2aa blocking peptide (CDBP5751)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	mfsd2aa major facilitator superfamily domain containing 2aa [Danio rerio (zebrafish)]
Official Symbol	mfsd2aa
Synonyms	mfsd2aa; major facilitator superfamily domain containing 2aa; NLS1-A; mfsd2a; zgc:101615; sodium-dependent lysophosphatidylcholine symporter 1-A; sodium-dependent LPC symporter 1-A; major facilitator superfamily domain containing 2a; major facilitator superfamily domain-containing protein 2A-A
Entrez Gene ID	492810
mRNA Refseq	NM_001007451
Protein Refseq	NP_001007452
UniProt ID	Q5U3U7

Function

phospholipid transporter activity; symporter activity
