



MFSD2B blocking peptide (CDBP5752)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MFSD2B major facilitator superfamily domain containing 2B [Homo sapiens (human)]
Official Symbol	MFSD2B
Synonyms	MFSD2B; major facilitator superfamily domain containing 2B; major facilitator superfamily domain-containing protein 2B
Entrez Gene ID	388931
mRNA Refseq	NM_001080473
Protein Refseq	NP_001073942
UniProt ID	A6NFX1