



MLIP blocking peptide (CDBP5757)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MLIP muscular LMNA-interacting protein [Homo sapiens (human)]
Official Symbol	MLIP
Synonyms	MLIP; muscular LMNA-interacting protein; C6orf142; muscle-enriched A-type lamin interacting protein; muscular-enriched A-type laminin-interacting protein
Entrez Gene ID	90523
mRNA Refseq	NM_001281746
Protein Refseq	NP_001268675
UniProt ID	Q5VWP3