



# MOSPD1 blocking peptide (CDBP5760)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |                                                            |
|---------------|------------------------------------------------------------|
| Conjugate     | Unconjugated                                               |
| Applications  | Used as a blocking peptide in immunoblotting applications. |
| Format        | Liquid                                                     |
| Concentration | 200 µg/mL                                                  |
| Size          | 0.05 mg                                                    |
| Preservative  | None                                                       |
| Storage       | -20°C                                                      |

## GENE INFORMATION

|                 |                                                                                             |
|-----------------|---------------------------------------------------------------------------------------------|
| Gene Name       | <a href="#">MOSPD1 motile sperm domain containing 1 [ Homo sapiens (human) ]</a>            |
| Official Symbol | MOSPD1                                                                                      |
| Synonyms        | MOSPD1; motile sperm domain containing 1; DJ473B4; motile sperm domain-containing protein 1 |
| Entrez Gene ID  | <a href="#">56180</a>                                                                       |
| mRNA Refseq     | <a href="#">NM_019556</a>                                                                   |
| Protein Refseq  | <a href="#">NP_062456</a>                                                                   |
| UniProt ID      | Q9UJG1                                                                                      |
| Function        | structural molecule activity                                                                |