



MTERF3 blocking peptide (CDBP5768)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MTERF3 mitochondrial transcription termination factor 3 [Homo sapiens (human)]
Official Symbol	MTERF3
Synonyms	MTERF3; mitochondrial transcription termination factor 3; CGI-12; MTERFD1; mTERF domain-containing protein 1, mitochondrial; MTERF domain containing 1
Entrez Gene ID	51001
mRNA Refseq	NM_001286643
Protein Refseq	NP_001273572
UniProt ID	Q96E29
Pathway	Mitochondrial Gene Expression

Function	transcription regulatory region DNA binding
-----------------	---
