



MTERF4 blocking peptide (CDBP5769)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MTERF4 mitochondrial transcription termination factor 4 [Homo sapiens (human)]
Official Symbol	MTERF4
Synonyms	MTERF4; mitochondrial transcription termination factor 4; MTERFD2; mTERF domain-containing protein 2; MTERF domain containing 2
Entrez Gene ID	130916
mRNA Refseq	NM_182501
Protein Refseq	NP_872307
UniProt ID	Q7Z6M4
Function	protein binding; rRNA binding