



MTERF2 blocking peptide (CDBP5770)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MTERF2 mitochondrial transcription termination factor 2 [Homo sapiens (human)]
Official Symbol	MTERF2
Synonyms	MTERF2; mitochondrial transcription termination factor 2; mTERFL; MTERFD3; mTERF domain-containing protein 3, mitochondrial; MTERF domain containing 3; mitochondrial transcription termination factor-like protein
Entrez Gene ID	80298
mRNA Refseq	NM_001033050
Protein Refseq	NP_001028222
UniProt ID	Q49AM1

Function

protein homodimerization activity; transcription regulatory region DNA binding
