



MYCT1 blocking peptide (CDBP5774)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MYCT1 myc target 1 [Homo sapiens (human)]
Official Symbol	MYCT1
Synonyms	MYCT1; myc target 1; MTLC; myc target protein 1; myc target in myeloid cells 1; myc target in myeloid cells protein 1
Entrez Gene ID	80177
mRNA Refseq	NM_025107
Protein Refseq	NP_079383
UniProt ID	Q8N699
Pathway	Validated targets of C-MYC transcriptional activation