



MYT1L blocking peptide (CDBP5778)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MYT1L myelin transcription factor 1-like [Homo sapiens (human)]
Official Symbol	MYT1L
Synonyms	MYT1L; myelin transcription factor 1-like; NZF1; ZC2HC4B; myelin transcription factor 1-like protein; myT1-L; neural zinc finger transcription factor 1
Entrez Gene ID	23040
mRNA Refseq	NM_015025
Protein Refseq	NP_055840
UniProt ID	Q9UL68
Function	DNA binding; sequence-specific DNA binding transcription factor activity; zinc ion binding