



NKIRAS1 blocking peptide (CDBP5780)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	NKIRAS1 NFKB inhibitor interacting Ras-like 1 [Homo sapiens (human)]
Official Symbol	NKIRAS1
Synonyms	NKIRAS1; NFKB inhibitor interacting Ras-like 1; KBRAS1; kappaB-Ras1; NF-kappa-B inhibitor-interacting Ras-like protein 1; kappa B-ras 1; kappa B-Ras protein 1; I-kappa-B-interacting Ras-like protein 1; NFKB inhibitor interacting Ras-like protein 1
Entrez Gene ID	28512
mRNA Refseq	NM_020345
Protein Refseq	NP_065078
UniProt ID	Q9NYS0

Pathway	TNF-alpha/NF-kB Signaling Pathway
Function	GTP binding; GTPase activity